

Apc Rabbit Anti-Human/Monkey Cd19 Mab - Cd19

Item Number: CM0010420

Introduction

APC-conjugated rabbit monoclonal antibody against CD19, designed for flow cytometry (FC) applications. Demonstrates cross-reactivity with human, cynomolgus, and rhesus monkey. Ideal for immunophenotyping and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID (Monkey)	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19 (XP_005591597.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus, Rhesus
Protein Weight	63kDa
Shipping	Ice bag