


ApC Rabbit Anti-Human/Monkey Cd197/Ccr7 Monoclonal Antibody For Fc - P32248

Item Number: CM0005225



Introduction

APC-conjugated rabbit monoclonal antibody targeting CD197/CCR7 (P32248), validated for flow cytometry. Recognizes human and cynomolgus monkey. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD197/CCR7 mAb
Remarks/Alias	BLR2; EBI1; CCR-7; CD197; CDw197; CMKBR7; CC-CCR-7
Species	Human
GeneID (Human)	1236
GeneID	1236
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P32248
Protein Weight	43kDa
Shipping	Ice bag