


## Apc Rabbit Anti-Human/Monkey Cd20 Mab - P11836

Item Number: CM0005338



### Introduction

APC-conjugated rabbit monoclonal antibody specific for human and cynomolgus CD20 (MS4A1, B1, Bp35, Leu-16). Validated for flow cytometry. Recognizes the B-lymphocyte antigen P11836, 33 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | APC Rabbit anti-Human/Monkey CD20 mAb                                                                |
| Remarks/Alias    | B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16                                                              |
| Species          | Human                                                                                                |
| GeneID (Human)   | 931                                                                                                  |
| GeneID           | 931                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 (NP_068769.2) |
| Source           | Rabbit                                                                                               |
| Category         | Monoclonal Antibodies                                                                                |
| Application      | FC                                                                                                   |
| Cross Reactivity | Human, Cynomolgus                                                                                    |
| SWISS            | P11836                                                                                               |
| Protein Weight   | 14kDa/33kDa                                                                                          |
| Shipping         | Ice bag                                                                                              |