

APC Rabbit Anti-Human/Monkey CD25 Monoclonal Antibody For Flow Cytometry - P01589

Item Number: CM0006233



Introduction

APC-conjugated rabbit anti-human/monkey CD25 monoclonal antibody for flow cytometry. Targets IL-2 receptor subunit alpha (p55, TAC antigen) on regulatory T cells; validated in human and cynomolgus.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
GeneID (Human)	3559
GeneID	3559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-213 of human CD25 (NP_000408.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag