


Apc Rabbit Anti-Human/Monkey Cd38 Polymab® For Flow Cytometry - P28907

Item Number: CM0012751



Introduction

APC-conjugated rabbit monoclonal antibody against human and cynomolgus CD38, validated for flow cytometry. Targets the CD38 ectoenzyme/receptor involved in calcium mobilization, lymphocyte activation, and cell adhesion. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD38 PolymAb®
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag