

ApC Rabbit Anti-Human/Monkey Cd69 Mab For Fc - Q07108

Item Number: CM0010465



Introduction

APC-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD69 (UniProt Q07108). Designed for flow cytometry, recognizes early activation antigen on activated T cells, ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID (Human)	969
GeneID	969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag