


ApC Rabbit Anti-Human/Monkey Cd86 Monoclonal Antibody For Fc - P42081

Item Number: CM0012479



Introduction

APC-conjugated rabbit monoclonal antibody targeting human/monkey CD86 (B7-2), validated for flow cytometry. Ideal for investigating T-cell costimulation and immune regulation in antigen-presenting cells.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD86 mAb
Remarks/Alias	B70; B7-2; B7.2; LAB72; CD28LG2
Species	Human
GeneID (human)	942
GeneID	942
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-241 of human CD86 (NP_787058.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P42081
Protein Weight	38kDa
Shipping	Ice bag