


Apc Rabbit Anti-Human/Monkey Cd8A Mab For Fc - P01732

Item Number: CM0008347



Introduction

APC-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD8a, validated for flow cytometry. Suitable for T-cell immunophenotyping and immune monitoring in diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8A; CD8; Leu2; MAL; p32; CD8a molecule
Species	Human
Gene ID (Human)	925
Gene ID	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag