

ApC Rabbit Anti-Human/Monkey Ki67 Monoclonal Antibody For Intracellular Flow Cytometry - P46013

Item Number: CM0012705



Introduction

APC-conjugated rabbit monoclonal antibody against human and cynomolgus Ki67 (P46013) for intracellular flow cytometry. Reliable proliferation marker for cell cycle analysis in oncology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey Ki67 mAb
Remarks/Alias	KIA; MIB-; MIB-1; PPP1R105
Species	Human
Gene ID (Human)	4288
Gene ID	4288
Immunogen	Synthetic peptide corresponding to amino acids 2271-2370 of human Ki67 (NP_002408.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human, Cynomolgus
SWISS	P46013
Protein Weight	359kDa
Shipping	Ice bag