

ApC Rabbit Anti-Human/Monkey Pd-1/Cd279 Mab For Flow Cytometry - Q15116

Item Number: CM0010216



Introduction

APC-conjugated rabbit monoclonal antibody targeting human and cynomolgus PD-1/CD279 (aa 25-168). Validated for flow cytometry, this antibody supports immuno-oncology research and immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human/Monkey PD-1/CD279 mAb
Remarks/Alias	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1
Species	Human
Gene ID (human)	5133
Gene ID	5133
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-168 of human PD-1/CD279 (NP_005009.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q15116
Protein Weight	32kDa
Shipping	Ice bag