

## ApC Rabbit Anti-Human/Mouse Pdgfra/Cd140A Monoclonal Antibody For Fc - P16234

Item Number: CM0004491



### Introduction

APC-conjugated rabbit monoclonal antibody targeting human and mouse PDGFR $\alpha$ /CD140a (SWISS: P16234, 123kDa). Validated for flow cytometry, ideal for cell proliferation, development, and oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                                              |
|------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | APC Rabbit anti-Human/Mouse PDGFR $\alpha$ /CD140a mAb                                                                             |
| Remarks/Alias    | CD140A; PDGFR2; PDGFR-2                                                                                                            |
| Species          | Human                                                                                                                              |
| GeneID (Human)   | 5156                                                                                                                               |
| GeneID           | 5156                                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human PDGFR $\alpha$ /CD140a (NP_006197.1) |
| Source           | Rabbit                                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                                              |
| Application      | FC                                                                                                                                 |
| Cross Reactivity | Human, Mouse                                                                                                                       |
| SWISS            | P16234                                                                                                                             |
| Protein Weight   | 123kDa                                                                                                                             |
| Shipping         | Ice bag                                                                                                                            |