

Apc Rabbit Anti-Mouse Cd103 Mab For Fc - Q60677

Item Number: CM0010402

Introduction

The APC Rabbit anti-Mouse CD103 mAb (clone aM290) is a monoclonal antibody conjugated to APC, validated for flow cytometry. It specifically recognizes mouse CD103 (Integrin alpha-E), a receptor for E-cadherin involved in intraepithelial lymphocyte homing. Ideal for mucosal immunology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD103 mAb
Remarks/Alias	CD103; aM290; alpha-E1; A530055J10; alpha-M290
Species	Mouse
GeneID (human)	3682
GeneID	16407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 191-385 of mouse CD103 (NP_032425.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60677
Protein Weight	129kDa
Shipping	Ice bag