

Apc Rabbit Anti-Mouse Cd124/Il-4Ra Monoclonal Antibody For Flow Cytometry - P16382

Item Number: CM0005486

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD124/IL-4Ra. Validated for flow cytometry to study IL-4/IL-13 signaling, Th2 differentiation, and allergic inflammation. - P16382

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD124/IL-4Ra mAb
Remarks/Alias	Il4r; CD124
Species	Mouse
GeneID (Human)	3566
GeneID	16190
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-233 of mouse CD124/IL-4Ra (NP_001008700.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16382
Protein Weight	87kDa
Shipping	Ice bag