

Apc Rabbit Anti-Mouse Cd138/Syndecan-1 Mab For Fc - P18828

Item Number: CM0011189

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD138/Syndecan-1 for flow cytometry. Recombinant immunogen (aa23-255). Specific for mouse samples. Suitable for membrane proteoglycan detection in exosome and matrix interaction studies.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Mouse CD138/Syndecan-1 mAb
Remarks/Alias	Sstn; Synd; CD138; Synd1; syn-1; CD138/Syndecan-1
Species	Mouse
GeneID (human)	6382
GeneID	20969
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-255 of mouse CD138/Syndecan-1 (NP_035649.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18828
Protein Weight	33kDa
Shipping	Ice bag