

Apc Rabbit Anti-Mouse Cd150/Slam Monoclonal Antibody For Fc - Q9Qum4

Item Number: CM0011442

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD150/SLAM (SLAM family member 1), validated for flow cytometry. Ideal for studying immune cell activation, differentiation, and SLAMF1 signaling in innate and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD150/SLAM mAb
Remarks/Alias	Slam; CD150; IPO-3; CDw150; ESTM51; 4933415F16
Species	Mouse
Gene ID (Human)	6504
Gene ID	27218
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-242 of mouse CD150/SLAM (NP_038758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QUM4
Protein Weight	36kDa/38kDa
Shipping	Ice bag