

Apc Rabbit Anti-Mouse Cd155/Pvr Mab - Q91Wp1

Item Number: CM0012448

Introduction

APC-conjugated rabbit anti-mouse CD155/Pvr monoclonal antibody for flow cytometry. Specific to mouse CD155, ideal for immune receptor research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD155/Pvr mAb
Remarks/Alias	PVS; mE4; HVED; Taa1; CD155; Tage4; necl-5; D7Ert458e; 3830421F03Rik
Species	Mouse
Gene ID (Human)	5817
Gene ID	52118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-348 of mouse CD155/Pvr (NP_612564.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q91WP1
Protein Weight	45kDa
Shipping	Ice bag