

Apc Rabbit Anti-Mouse Cd21/Cd35 Monoclonal Antibody For Fc - P19070

Item Number: CM0007899

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD21/CD35 (CR2/C3DR). Validated for flow cytometry, reacts with mouse B-cells. Useful for complement and B-cell immunology studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD21/CD35 mAb
Remarks/Alias	Cr1; C3DR; CD21; CD35; Cr-1; Cr-2
Species	Mouse
GenelD (human)	1380
GenelD	12902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-274 of mouse CD21/CD35 (NP_031784.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P19070
Protein Weight	113kDa
Shipping	Ice bag