

Apc Rabbit Anti-Mouse Cd22 Monoclonal Antibody For Fc - P35329

Item Number: CM0009913

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD22 (Siglec-2, BL-CAM) for flow cytometry. Validated for mouse reactivity, ideal for B-cell immunophenotyping and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
GeneID (human)	933
GeneID	12483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa/93kDa/97kDa
Shipping	Ice bag