

Apc Rabbit Anti-Mouse Cd24 Monoclonal Antibody For Flow Cytometry - P24807

Item Number: CM0005930

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD24 (HSA), validated for flow cytometry. Suitable for immune cell differentiation, B-cell activation, autoimmunity, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
Gene ID (Human)	100133941
Gene ID	12484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag