

APC Rabbit Anti-Mouse Cd272/Btla Monoclonal Antibody For Flow Cytometry - Q7Tsa3

Item Number: CM0004665



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD272/BTLA, an inhibitory receptor on lymphocytes. Designed for flow cytometry (FC) to study immune regulation in research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD272/BTLA mAb
Remarks/Alias	A630002H24
Species	Mouse
Gene ID (Human)	151888
Gene ID (Mouse)	208154
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q7TSA3
Protein Weight	22kDa/34kDa
Shipping	Ice bag