


Apc Rabbit Anti-Mouse Cd304/Neuropilin-1 Monoclonal Antibody For Flow Cytometry - P97333

Item Number: CM0011571



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse Neuropilin-1 (CD304), validated for flow cytometry. Involved in angiogenesis, neuronal development, and organogenesis. Ideal for cardiovascular and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD304/Neuropilin-1 mAb
Remarks/Alias	Nrp; NP-1; Npn1; NPN-1; C530029I03
Species	Mouse
GeneID (Human)	8829
GeneID (Mouse)	18186
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-856 of mouse CD304/Neuropilin-1 (NP_032763.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P97333
Protein Weight	103kDa
Shipping	Ice bag