


Apc Rabbit Anti-Mouse Cd317/Bst2 Monoclonal Antibody For Fc - Q8R2Q8

Item Number: CM0007731



Introduction

APC-conjugated rabbit anti-mouse CD317/BST2 monoclonal antibody for flow cytometry. Recognizes mouse bone marrow stromal antigen 2 (19kDa). Useful for antiviral restriction studies, innate immunity, and plasmacytoid dendritic cell research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD317/BST2 mAb
Remarks/Alias	GREG; Bst-2; CD317; DAMP-1; 2310015I10Rik
Species	Mouse
Gene ID (Human)	684
Gene ID	69550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317/BST2 (NP_932763.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8R2Q8
Protein Weight	19kDa
Shipping	Ice bag