

Apc Rabbit Anti-Mouse Cd32 Monoclonal Antibody For Flow Cytometry (Fc) - P08101

Item Number: CM0009484

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD32 (FcγRIIB). Specific for flow cytometry (FC) analysis of low affinity IgG Fc receptor involved in immune regulation. Suitable for immunology and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD32 mAb
Remarks/Alias	CD32; Fcgr2; Fcr-2; Fcr-3; Ly-17; LyM-1; Lym-1; fcRII; FcgRII; Fcgr2a; Ly-m20; Fc[g]RII; F630109E10Rik
Species	Mouse
GeneID (Human)	2213
GeneID	14130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of mouse CD32 (NP_034317.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08101
Protein Weight	31kDa/33kDa/36kDa
Shipping	Ice bag