

Apc Rabbit Anti-Mouse Cd34 Monoclonal Antibody For Flow Cytometry - Q64314

Item Number: CM0012071

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD34, specifically designed for flow cytometry detection. Targets hematopoietic progenitor cells with a calculated molecular weight of ~41 kDa. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Mouse
GeneID (Human)	947
GeneID	12490
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41kDa
Shipping	Ice bag