

APC Rabbit Anti-Mouse Cd366 (Tim-3) Monoclonal Antibody For Flow Cytometry - Q8Vim0

Item Number: CM0009975



Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD366 (TIM-3) for flow cytometry. Ideal for studying T-cell exhaustion, immune checkpoints, and macrophage activation. Validated in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD366 (TIM-3) mAb
Remarks/Alias	Tim3; TIM-3; Timd3
Species	Mouse
GeneID (Human)	84868
GeneID	171285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-191 of mouse CD366 (TIM-3) (NP_599011.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8VIM0
Protein Weight	22kDa/31kDa
Shipping	Ice bag