

Apc Rabbit Anti-Mouse Cd38 Monoclonal Antibody For Flow Cytometry - P56528

Item Number: CM0007885



Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD38 (ADPRC 1, I-19). Reacts with mouse, suitable for flow cytometry analysis of CD38 expression in immunology and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD38 mAb
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1
Species	Mouse
GeneID (human)	952
GeneID	12494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-304 of mouse CD38 (NP_031672.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag