

Apc Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Fc - P06332

Item Number: CM0007758

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD4 (L3T4, Ly-4), validated for flow cytometry (FC). Suitable for identifying T-helper cells and other CD4+ immune populations in murine samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID	12504
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag