

APC-Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Flow Cytometry - P15379

Item Number: CM0011998



Introduction

APC--conjugated rabbit anti-mouse CD44 monoclonal antibody for flow cytometry. Reacts with mouse CD44 (Ly-24, Pgp-1, HERMES), a hyaluronate receptor involved in lymphocyte homing, adhesion, and inflammation. Suitable for mouse research samples.

[Learn More](#)

Parameter	Value
Product Name	APC-Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	86kDa
Shipping	Ice bag