

Apc Rabbit Anti-Mouse Cd47 Mab - Q61735

Item Number: CM0010041

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD47 (IAP), validated for flow cytometry. Ideal for studying cell adhesion, immune regulation, and thrombospondin signaling. Suitable for research in immunology and cell biology.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD47 mAb
Remarks/Alias	IAP; Itgp; 9130415E20Rik; B430305P08Rik
Species	Mouse
GeneID (Human)	961
GeneID	16423
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of mouse CD47 (NP_034711.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61735
Protein Weight	33kDa
Shipping	Ice bag