

Apc Rabbit Anti-Mouse Cd59A Monoclonal Antibody For Fc - O55186

Item Number: CM0011503

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD59a, a potent MAC inhibitor. For flow cytometry applications in mouse samples. Suitable for complement, immunology, and therapeutic research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD59a mAb
Remarks/Alias	Cd59; MACIF; MAC-IP
Species	Mouse
GeneID	12509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-96 of mouse CD59a (NP_031678.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O55186
Protein Weight	14 kDa
Shipping	Ice bag