

Apc Rabbit Anti-Mouse Cd71 Monoclonal Antibody For Fc - Q62351

Item Number: CM0008611

Introduction

High-quality APC-conjugated rabbit monoclonal antibody against mouse CD71 (Transferrin receptor protein 1) validated for flow cytometry. Ideal for studying iron uptake and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD71 mAb
Remarks/Alias	TR; TFR; p90; CD71; TFR1; Trfr; Mtvr1; 2610028K12Rik; E430033M20Rik
Species	Mouse
GeneID (Human)	7037
GeneID	22042
Immunogen	Recombinant protein corresponding to a sequence within amino acids 89-763 of mouse CD71 (NP_035768.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62351
Protein Weight	86kDa
Shipping	Ice bag