

Apc Rabbit Anti-Mouse Cd86 Monoclonal Antibody For Flow Cytometry - P42082

Item Number: CM0009635

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD86, designed for flow cytometry (FC). Recognizes mouse CD86 (SWISS P42082), ~35kDa. Ideal for T-cell activation and immune regulation studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD86 mAb
Remarks/Alias	B7; B70; MB7; B7-2; B7.2; CLS1; Ly58; ETC-1; Ly-58; MB7-2; Cd28l2; TS/A-2; CD86
Species	Mouse
Gene ID (human)	942
Gene ID	12524
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-245 of mouse CD86 (NP_062261.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P42082
Protein Weight	35kDa
Shipping	Ice bag