

Apc Rabbit Anti-Mouse Cd86 Monoclonal Antibody For Flow Cytometry - P42082

Item Number: CM0008264

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD86 (B7-2), validated for flow cytometry. Targets amino acids 24-245 of mouse CD86. Ideal for immunology research on T-cell activation and costimulation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD86 mAb
Remarks/Alias	B7, B70, MB7, B7-2, B7.2, CLS1, Ly58, ETC-1, Ly-58, MB7-2, Cd28l2, TS/A-2, CD86
Species	Mouse
GenelD (Human)	942
GenelD	12524
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-245 of mouse CD86 (NP_062261.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P42082
Protein Weight	35 kDa
Shipping	Ice bag