

Apc Rabbit Anti-Mouse Cd93 Monoclonal Antibody For Flow Cytometry - O89103

Item Number: CM0004720

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD93 (C1qR, AA4). Suitable for flow cytometry (FC). Enables study of phagocytosis, endothelial function, and angiogenesis. Offers batch-to-batch reproducibility.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD93 mAb
Remarks/Alias	Ly68; AA4.1; C1qr1; C1qrp; 6030404G09Rik
Species	Mouse
Gene ID (Human)	22918
Gene ID	17064
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 23-181 of mouse CD93 (NP_034870.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O89103
Protein Weight	69kDa
Shipping	Ice bag