


ApC Rabbit Anti-Mouse Cd96/Tactile Monoclonal Antibody For Fc - Q3U0X8

Item Number: CM0011436



Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD96/TACTILE for flow cytometry. Validated for mouse samples. Targets T-cell surface protein involved in NK cell adhesion during immune response.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Mouse CD96/TACTILE mAb
Remarks/Alias	Tactile; 1700109I12Rik
Species	Mouse
GeneID (Human)	10225
GeneID	84544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-536 of mouse CD96/TACTILE (NP_115854.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q3U0X8
Protein Weight	67kDa
Shipping	Ice bag