

APC Rabbit Anti-Mouse H-2Kd/H-2Dd Monoclonal Antibody For Fc - P01902/P01900

Item Number: CM0012167



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse H-2Kd/H-2Dd MHC class I molecules. Validated for flow cytometry, ideal for immune cell typing, antigen presentation studies, and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse H-2Kd/H-2Dd mAb
Remarks/Alias	H2-D1, H-2D; H2-K1, H-2K, H2-K
Species	Mouse
Gene ID (Mouse)	14972/14964
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-301 of mouse H-2Kd/H-2Dd (NP_612564.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01902/P01900
Protein Weight	41kDa
Shipping	Ice bag