

Apc Rabbit Anti-Mouse Icos/Cd278 Monoclonal Antibody For Fc - Q9Wvs0

Item Number: CM0011084

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse ICOS/CD278 (UniProt Q9WVS0), validated for flow cytometry. Suitable for studying T-cell costimulation, follicular helper cell differentiation, and related immunological assays.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse ICOS/CD278 mAb
Remarks/Alias	H4; CCLP; AILIM; CRP-1; Ly115
Species	Mouse
GeneID (Human)	29851
GeneID	54167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-142 of mouse ICOS/CD278 (NP_059508.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WVS0
Protein Weight	23KDa
Shipping	Ice bag