

Apc Rabbit Anti-Mouse Igg1 Monoclonal Antibody For Fc - P01868

Item Number: CM0004892



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse IgG1 constant region. Designed for flow cytometry (FC) applications, this antibody specifically detects mouse IgG1 in research and diagnostic workflows. Ideal for immunoassay development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse IgG1 mAb
Remarks/Alias	IgG1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-324 of mouse IgG1(P01868)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01868
Protein Weight	36kDa
Shipping	Ice bag