

Apc Rabbit Anti-Mouse Igg2B Monoclonal Antibody For Fc - P01867-2

Item Number: CM0006496

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse IgG2b heavy chain constant region. Validated for flow cytometry (FC). Reacts with mouse. Suitable for immunology research, diagnostic assay development, and B cell analysis.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse IgG2b mAb
Remarks/Alias	IgG2b; Igh-3; gamma2b
Species	Mouse
Gene ID	16016
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 14-113 of mouse IgG2b (P01867-2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01867-2
Protein Weight	37kDa
Shipping	Ice bag