

Apc Rabbit Anti-Mouse Il-13 Monoclonal Antibody For Intracellular Flow Cytometry - P20109

Item Number: CM0011144

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse IL-13, validated for intracellular flow cytometry. Suitable for studying allergic inflammation, Th2 responses, and IL-13 signaling pathways in research models.

[Learn More](#)

Parameter	Specification
Product Name	APC Rabbit anti-Mouse IL-13 mAb
Remarks/Alias	IL-13
Species	Mouse
GeneID (Human)	3596
GeneID (Mouse)	16163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-131 of mouse IL-13 (NP_032381.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P20109
Protein Weight	14kDa
Shipping	Ice bag