

Apc Rabbit Anti-Mouse Il-17A Monoclonal Antibody For Flow Cytometry (Intracellular) - Q62386

Item Number: CM0005593

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse IL-17A, validated for intracellular flow cytometry. Reacts with mouse IL-17A, a key Th17 cytokine involved in inflammation, host defense, and epithelial barrier integrity. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse IL-17A mAb
Remarks/Alias	Il17; Ctla8; IL-17; Ctla-8; IL-17A
Species	Mouse
GeneID (Human)	3605
GeneID	16171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse IL17A. (NP_034682.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q62386
Protein Weight	17 kDa
Shipping	Ice bag