

Apc Rabbit Anti-Mouse IL1R1/Cd121A Monoclonal Antibody For Flow Cytometry - P13504

Item Number: CM0011941

Introduction

APC-conjugated rabbit monoclonal antibody against mouse IL1R1/CD121a, validated for flow cytometry. Ideal for profiling interleukin-1 receptor expression in immune and neuronal research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse IL1R1/CD121a mAb
Remarks/Alias	IL-iR; CD121a; CD121b; IL-1R1; Il1r-1; IL-1R-1; IL-1RT1; IL-1RT-1; IL-1R-alpha
Species	Mouse
GeneID (Human)	3554
GeneID	16177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-338 of mouse IL1R1/CD121a (NP_032388.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13504
Protein Weight	67kDa
Shipping	Ice bag