


ApC Rabbit Anti-Mouse Klrg1/Mafa Monoclonal Antibody For Fc - O88713

Item Number: CM0006006



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse KLRG1/MAFA, suitable for flow cytometry. Recombinant immunogen covering amino acids 57-188. Expressed on NK cells and activated CD8 T-cells; involved in immune regulation.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse KLRG1/MAFA mAb
Remarks/Alias	MAFA; 2F1-Ag; MAFA-L
Species	Mouse
GeneID (Human)	10219
GeneID	50928
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 57-188 of mouse KLRG1/MAFA (NP_058666.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88713
Protein Weight	21kDa
Shipping	Ice bag