

Apc Rabbit Anti-Mouse Ly-6C Monoclonal Antibody For Flow Cytometry - P0Cw02

Item Number: CM0009085

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse Ly-6C (Ly6c1). Validated for flow cytometry (FC), specific to mouse, with a ~15 kDa protein weight. Ideal for immune cell profiling research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse Ly-6C mAb
Remarks/Alias	Ly6c; Ly-6C; Ly-6C1
Species	Mouse
Gene ID	17067
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of Mouse Ly-6C (NP_001238984.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P0CW02
Protein Weight	15 kDa
Shipping	Ice bag