

ApC Rabbit Anti-Mouse Podoplanin Monoclonal Antibody For Flow Cytometry - Q62011

Item Number: CM0008280



Introduction

APC-conjugated rabbit anti-mouse Podoplanin monoclonal antibody for flow cytometry. Validated for mouse samples; immunogen covers amino acids 21-133. Ideal for lymphatic, tumor, and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse Podoplanin mAb
Remarks/Alias	E11; T1a; Gp38; OTS-8; T1alpha; RANDAM-2; T1-alpha
Species	Mouse
GeneID (human)	10630
GeneID (mouse)	14726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-133 of mouse Podoplanin (NP_034359.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62011
Protein Weight	18kDa
Shipping	Ice bag