


ApC Rabbit Anti-Mouse Tim-1/Kim-1/Havcr Monoclonal Antibody For Fc - Q5Qns5

Item Number: CM0005909



Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse TIM-1/KIM-1/HAVCR for flow cytometry. Recombinant immunogen (aa 22-237) ensures specificity. Suitable for immunology, kidney injury, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse TIM-1/KIM-1/HAVCR mAb
Remarks/Alias	Tim1; KIM-1; TIM-1; Timd1; TIM-1/KIM-1/HAVCR
Species	Mouse
GeneID (human)	26762
GeneID	171283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-237 of mouse TIM-1/KIM-1/HAVCR (NP_599009.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q5QNS5
Protein Weight	33kDa
Shipping	Ice bag