

ApC Rabbit Anti-Mouse Tnf-A Monoclonal Antibody For Fc (Intra) - P06804

Item Number: CM0011495



Introduction

APC-conjugated rabbit monoclonal antibody specific for mouse TNF- α , validated for intracellular flow cytometry. Detects mouse TNF- α (cachectin, TNFSF2) with high affinity. Ideal for inflammation, immunology, and cytokine research applications.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse TNF- α mAb
Remarks/Alias	DIF; Tnfa; TNF-a; TNFSF2; Tnlg1f; Tnfsf1a; TNFalpha; TNF-alpha
Species	Human
GeneID (Human)	7124
GeneID	21926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-235 of mouse TNF- α (NP_038721.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P06804
Protein Weight	26kDa
Shipping	Ice bag