

Apc Rabbit Anti-Mouse/Human Cd55 Monoclonal Antibody For Fc - Q61475

Item Number: CM0008724

Introduction

APC-conjugated rabbit monoclonal antibody recognizing mouse and human CD55 (decay-accelerating factor), a GPI-anchored complement inhibitor. Suitable for flow cytometry-based analysis of complement regulatory pathways. UniProt Q61475.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse/Human CD55 mAb
Remarks/Alias	Daf; Daf1; Daf-GPI; GPI-DAF
Species	Mouse
GeneID (Human)	1604
GeneID	13136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-361 of mouse CD55/DAF (NP_034146.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q61475
Protein Weight	42KDa
Shipping	Ice bag