

Apc Rabbit Anti-Mouse/Rat Cd54 (Icam-1) Monoclonal Antibody For Flow Cytometry - P13597

Item Number: CM0009747

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse/rat CD54 (ICAM-1) for flow cytometry. Validated for mouse and rat cross-reactivity. Ideal for leukocyte adhesion and migration research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse/Rat CD54 mAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2
Species	Mouse
GeneID (Human)	3383
GeneID	15894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse CD54 (NP_034623.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P13597
Protein Weight	59kDa
Shipping	Ice bag