

Apc Rabbit Igg Isotype Control - Apc Rabbit Igg Isotype Control

Item Number: CM0008133

Introduction

APC-conjugated Rabbit IgG isotype control monoclonal antibody for flow cytometry (FC) applications. Ideal as a matched negative control in immunofluorescence staining panels, ensuring assay specificity and background assessment.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit IgG isotype control
Immunogen	Synthetic peptide (A synthesized peptide derived from rabbit IgG isotype control)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag